

# RISIO<sup>®</sup>

whitening



**Clinical protocol  
for effective  
deep bleaching**

## What is deep bleaching?

Using a combination of Hydrogen and Carbamide peroxides over a sustained period. The tooth is prepared for whitening using HP (Hydrogen peroxide). This works by revitalizing the tooth's ability to absorb oxygen and be extensively bleached by CP (Carbamide peroxide).

### 6% HP


Hydrogen peroxide breaks down into water and peroxide ions. These peroxide ions then bind to the chromophores which breaks the double bond. The single bonded chromophores absorb less light giving a whiter appearance preparing the tooth for treatment and familiarising the tooth with the process.

### 16% CP

Carbamide peroxide breaks down into hydrogen peroxide, urea and ammonia. The longer wear time allows the ingredient to penetrate the tooth microstructure deeper than that of HP resulting in better longer lasting results. Combination whitening has proven to be more effective over time than individual 1 hour or overnight whitening completed independently.


## Extrinsic and Intrinsic Stains

Colour molecules called chromophores bond together causing the tooth to look darker and darker. The colour molecules collect on the surface of the teeth (extrinsic stains) and the inside of the teeth (intrinsic stains).




## Primary and Secondary Factors of Teeth Whitening

The primary factors of teeth whitening are the concentration of the peroxide and the amount of time it is left on the teeth for absorption and breaking down the stain. The secondary factors increase the efficacy of the whitening oxygen with pH, temperature increase and preventing the escape of the whitening oxygen within a closed environment.



## Strength

Using the strongest available concentrations, we are ensuring that not only is your treatment as effective as possible but results will be long lasting. Historically, denser concentration coincides with more adverse effects, however, we take all necessary measures to ensure that adverse effects are minimised by the products of this treatment.



## Neutral pH

Our gels maintain a pH level of around 6 or slightly above. This stability ensures a highly effective yet comfortable treatment. Most companies have a whitening gel pH neutral on manufacture but quickly degrades in transit. Our gels are continually pH neutral and are regularly examined for efficacy. We only want the gel activating on the tooth surface.

## Treatment Protocol

1

Preparing the teeth is vital for a comfortable treatment. The patient should use the sensitivity relief serum before whitening. This will protect the teeth from adverse effects without affecting the treatment outcome. This can be used on demand.

2

On day one, apply the 6% HP application for one hour. This will prepare the teeth for longer and deeper applications of 16% CP.

3

Teeth are now prepared and ready for a longer denser treatment. Wearing the 16% CP overnight for 13 consecutive days results in substantial shade improvement.

4

On final day of treatment, the 6% HP is used for one hour. This finalises results and ensures an ultimate shade increase. There is assured treatment left for top ups going forward. At this point, we should be at shade B1. If not, we will provide additional gel free of charge.





Using this deep bleaching protocol alongside our relief serum, your patients are guaranteed to dramatically improve tooth shade without adverse sensitivity.





# RISIO<sup>®</sup>

**Contact us on:**

**E:** [info@risio.com](mailto:info@risio.com)

**T:** +44 (0) 191 8237990

**W:** [risio.com](http://risio.com)

ALS Dental Registered Office:  
85 Great Portland Street,  
London, England, W1W 7LT

**Find out more**

