

## CLINICAL PARAMETERS

1

1mm maximum  
centre line  
movement

2

2mm maximum  
buccal expansion

3

1.5mm maximum  
intrusion and  
extrusion

4

5-5 movement only

5

Up to 30-degree  
rotation

6

Over correction  
of the last 3 aligners



W: [als-dental.com](http://als-dental.com)  
E: [risio@als-dental.com](mailto:risio@als-dental.com)

Ashford Orthodontics Ltd. Unit 14B Southwick  
Industrial Estate Riverside Road, Sunderland SR5 3JG.

FIND OUT MORE

